

Defect Review in the Photonics Revolution

Aaron Lewis

Nanonics Imaging Ltd.

The Manhat Technology Park

Malcha, Jerusalem 91487 ISRAEL


Tel: 972-2-678-9573

Fax: 972-2-648-0827

Email: info@nanonics.co.il Website: www.nanonics.co.il

Report Documentation Page				Form Approved OMB No. 0704-0188	
Public reporting burden for the collection of information is estimated to average 1 hour per response, including the time for reviewing instructions, searching existing data sources, gathering and maintaining the data needed, and completing and reviewing the collection of information. Send comments regarding this burden estimate or any other aspect of this collection of information, including suggestions for reducing this burden, to Washington Headquarters Services, Directorate for Information Operations and Reports, 1215 Jefferson Davis Highway, Suite 1204, Arlington VA 22202-4302. Respondents should be aware that notwithstanding any other provision of law, no person shall be subject to a penalty for failing to comply with a collection of information if it does not display a currently valid OMB control number.					
1. REPORT DATE 18 APR 2000		2. REPORT TYPE N/A		3. DATES COVERED -	
4. TITLE AND SUBTITLE Defect Review in the Photonics Revolution				5a. CONTRACT NUMBER	
				5b. GRANT NUMBER	
				5c. PROGRAM ELEMENT NUMBER	
6. AUTHOR(S)				5d. PROJECT NUMBER	
				5e. TASK NUMBER	
				5f. WORK UNIT NUMBER	
7. PERFORMING ORGANIZATION NAME(S) AND ADDRESS(ES) Nanonics Imaging Limited				8. PERFORMING ORGANIZATION REPORT NUMBER	
9. SPONSORING/MONITORING AGENCY NAME(S) AND ADDRESS(ES)				10. SPONSOR/MONITOR'S ACRONYM(S)	
				11. SPONSOR/MONITOR'S REPORT NUMBER(S)	
12. DISTRIBUTION/AVAILABILITY STATEMENT Approved for public release, distribution unlimited					
13. SUPPLEMENTARY NOTES DARPA/MTO, WDM for Military Platforms Workshop held in McLean, VA on April 18-19, 2000, The original document contains color images.					
14. ABSTRACT					
15. SUBJECT TERMS					
16. SECURITY CLASSIFICATION OF:			17. LIMITATION OF ABSTRACT UU	18. NUMBER OF PAGES 20	19a. NAME OF RESPONSIBLE PERSON
a. REPORT unclassified	b. ABSTRACT unclassified	c. THIS PAGE unclassified			

Components of The DWDM Revolution



A New Direction for WDM Test & Measurement

- The Need - Decreasing Size and Increasing Integration
- The Technology - Near-field Optics
- 3. The level of Optical Resolution - 0.05 microns
- 4. New Correlations -
Light distribution

with simultaneous nanometer information on

Topography

Polarization


Electrical Characterization

Thermal Characterization

What is Near-Field Optics?


1. Collection or illumination of light through a sub-wavelength aperture

2. Scanning of sample or aperture relative to one another in the near-field




Near-field Optical Element

1. Tapered cantilevered metal-coated optical fiber probe
2. Simultaneous optical and topographical imaging
3. Simultaneous electrical [resistance, capacitance etc.] imaging



Deflection Force Sensing



The standard optical method

1. Cantilevered optical fibers are excellent atomic force sensors that can operate to detect surface topography either in contact or by being close but not in physical contact with the surface
2. Nanonics also provides non-optical means of deflection sensing

Optical Spectrum Analysis
Time Domain Measurements
Electrical Measurements
Polarization Analysis
Thermal Measurements

**CCD
CAMERA**

LENS

TIP SCANNER

**CANTILEVERED
OPTICAL FIBER**

**3D FLAT
SCANNER™**


SAMPLE SCANNER

LENS


**CCD
CAMERA**

Nanonics
Near-field/Far-field
Defect Review Station

System Flexibility



Placement of waveguide
for high resolution
injection of light at the
edge of the guide



Placement of the waveguide
for imaging evanescent
fields




Near-Field Optics (NSOM) Plays A Bridging Role Between Conventional Optical Microscopy And Atomic Force Microscopy



Demonstrating the uses of near-field optics


- Slab waveguide lasers
- V-Groove lasers
 - Correlating structure with light emission
 - Correlating mode structure with wavelength
 - Correlating mode structure with heat
- Optical waveguides
 - Optical fiber mode distribution with alterations in coupling
 - Polarization dependence of evanescent fields
 - Imaging the Tien effect
 - Imaging star couplers

Light distribution analysis with 0.05 micron spatial resolution of slab waveguide lasers emitting at 1.5 microns




A comparison of the (a) far-field and (b) near-field light distribution







***NSOM Light
Distribution
From The Laser
Cavity With An
Injected Current
That Is Below
The Threshold
For Lasing
Action***



***NSOM Light
Distribution
From The Laser
Cavity With An
Injected Current
That Is Above
The Threshold
For Lasing
Action***






***Injected Charge
Distribution
Measured Using
the Simultaneous
Atomic Force
Capabilities With
an Injected
Current Above
the Threshold for
Lasing Action***




*An Electron
Micrograph
Of The V
Groove Laser
Structure*

*Diagrammatic
Representation Of
The Structure
(B).*




1. Correlation of the light distribution and geometric structure of the v groove laser
2. Notice the 150 nm offset




Correlation of mode structure with wavelength alteration

QWR293 -NSOM








QWR294 -Temperature




Correlation of light distribution with thermal characteristics

Optical Fiber Output Analysis




Correlation of light distribution with fiber coupling




10 nm deep
topographic
alteration correlated
with light leakage
from a waveguide
that corresponds to
an ~ 0.05 dB guided
power loss

*Images taken
from Applied
Physics
Letters
Vol. 73,
1035-1037
(1998)*



Near-field
optical images of
the evanescent
field for TE and
TM polarization
of a
semiconductor
waveguide



Near-field optical image of the
star coupler section of a phasor
device

Summary

Near-field Optics

exceptional integrated information on today's components
with a resolution and an information content critical
for components in sight for tomorrow
and

being considered for tomorrow

[S. Shanhui, I. Appelbaum & J. D. Joannopoulos, "Near-field scanning optical microscopy as a simultaneous probe of fields and band structure of photonic crystals," Applied Physics Letters **75**, 3461 (1999)]

Nanonics Imaging Ltd.

The Manhat Technology Park
Malcha, Jerusalem 91487 ISRAEL

Tel: 972-2-678-9573

Fax: 972-2-648-0827

Email: info@nanonics.co.il Website: www.nanonics.co.il